



25 colours are available for your **PLAY SAFE** type.

Please indicate colour numbers:



Playsafe mouthguard out of up to **4** colours!
You can compose your own multi-coloured mouthguard out of the colours **1** to **15**.

Please indicate the colour numbers:

Example: 4 8 9 4



Please check the requested Playsafe type:

PLAY SAFE triple + **PLAY SAFE triple light**



Playsafe triple and triple light are suited for **any kind of sport**. Three-layered sports mouthguard with hard intermediate layer. Triple light is suited in case of cramped oral situation.



PLAY SAFE heavy-pro, 3-schichtig, with hard intermediate layer, for: **Field Hockey • Rollerblading • Squash • Polo • Ice Hockey • Karate • Kickboxing • Rugby** etc. In case of hard and pointed impacts.



PLAY SAFE light-pro, 3 layers, with hard intermediate layer. Suited in case of cramped oral situation. Same application area as Playsafe heavy-pro.



PLAY SAFE medium, 2 layers, for: **Boxing • Biking • Basketball • Soccer • Judo • Wrestling** etc. Sports where objects with a large surface area may be received by the jaw.



PLAY SAFE light, 2 layers, in case of cramped oral situation, for: **see medium**. Playsafe light protects opposing teeth hitting against each other from intraoral injuries.

Dear Athlete,
Every Playsafe mouthguard is fabricated for you by a trained person according to the Playsafe guidelines.

- ▶ You will get your Playsafe mouthguard in your dental practice. There an impression of your upper and lower jaw will be taken.
- ▶ In the dental lab, the plaster model and the sports mouthguard will then be produced.
- ▶ The control of fit and handing over will again take place in the dental practice.

An original Playsafe sports mouthguard is always laminated (min. two layers) and shows the Playsafe type and the name of the athlete.

Your Playsafe team:

Playsafe 4u: free individualisation of a Playsafe



Procedure:

- Choose preferred design and determine where it should be applied. Mostly at the visible surface during wearing. Because of the bite changing it is not recommended to apply designs at the occlusal surface.
- Erkodent needs the data of the design as .jpg or .pdf. The print will be effected immediately and sent by letter to the client.

All about Playsafe 4u auf www.erkodent.com

ERKODENT

Erkodent Erich Kopp GmbH • Siemensstraße 3 • 72285 Pfalzgrafenweiler Germany • Tel. + 49 (0) 74 45/85 01-0 • Fax + 49 (0) 74 45/20 92
info@erkodent.com • www.erkodent.com
ISO 9001:2008 • ISO 13485:2012

PLAY SAFE



Sports mouthguard

PLAY SAFE +
PLAY SAFE triple

ERKODENT

EN